

# IC15 ACTIVEX™ MATTRESS


Feel Rating Firm

Mattress Height 150mm

Risk Category High Risk

User Weight 45kg - 300kg

Warranty 10 Year Warranty



# **BED CODES & SIZES**

LONG SINGLE Code: IC15LS

Mattress Size: 2030 x 900

KING SINGLE Code: IC15KS

Mattress Size: 2030 x 1070

LONG DOUBLE Code: IC15LD

Mattress Size: 2030 x 1350

QUEEN Code: IC15Q

Mattress Size: 2030 x 1520

SPLIT QUEEN Code: IC15SQ

Mattress Size: 2030 x 750 (per side)

# **KEY FEATURES**

- 150mm Mattress Height
- 300kg Weight Rating
- · Zip off washable covers
- · Anti Slip Base Layer
- Contour Latch System
- Mutli layer immersion
- · Compatible with 15cm Zerotec waterproof covers
- Flip free design
- · Designed for use on adjustable beds
- 20+ waterlow rating
- 10 Year Warranty

#### **SPECIFICATIONS**

- ActiveX<sup>TM</sup> pressure relieving material
- · Vant ventilated sides and base fabric
- Bluetec temperature regulating cover
- · Built in firm edge support
- Flip free design
- · Zip off machine washable covers
- 45 300kg user weight
- Standard cover MVTR: NA.
- Medical Cover MVTR (g/24h·m2): 1635.0 Average
- Medical Cover BS7175: Crib 5
- ARTG Identifier: 259893

## PRESSURE CARE

Risk Category High

This rating is based on our average peak pressure of <24mmHg. For some high risk clients lcare Zerotec covers are recommended.

Icare risk category ratings are based on average peak pressure, Therapist recommendation, case studies, industry comparison and the Waterlow scale. For some high risk clients Icare Zerotec covers are recommended. The best way to assess a client on a pressure reducing surface is to monitor pressure points with visual checks over a period of 5 days.

## **CARE & CLEANING**

To launder the cover in a washing machine, set the main wash to a minimum of 65°C for no less than 10 minutes. The cover can be tumble dried on low heat and the temperature must not exceed 60°C. The foam can be autoclaved at 134°C.

