


# M2 M-LINE™ MATTRESS

Feel Rating	Medium
Mattress Height	150mm
Risk Category	High Risk
User Weight	40kg - 215kg
Warranty	5 Year Warranty



## BED CODES & SIZES

### LONG SINGLE

Code: ICM2LS  
Mattress Size: 2030 x 900

### KING SINGLE

Code: ICM2KS  
Mattress Size: 2030 x 1070

### LONG DOUBLE

Code: ICM2LD  
Mattress Size: 2030 x 1350

### QUEEN

Code: ICM2Q  
Mattress Size: 2030 x 1520

### SPLIT QUEEN

Code: ICM2SQ  
Mattress Size: 2030 x 760 (per side)

## KEY FEATURES

- 150mm Mattress Height
- 215kg Weight Rating
- Reduce roll-out
- Anti bottom out layer
- Waterproof Easy clean
- FR rated cover
- Flip free design
- Designed for use on adjustable beds
- 20+ waterlow rating
- 5 Year Warranty

## SPECIFICATIONS

- XP Support base and sides
- ZeroTec waterproof cover
- Built in firm edge support
- Flip free design
- Zip off machine washable 4 way stretch covers (replacement covers available if required)
- 40 - 215kg user weight
- Medical Cover MVTR (g/24h·m2): 1635.0 Average
- Medical Cover BS7175: Crib 5
- ARTG Identifier: 375851

## PRESSURE CARE

Risk Category High

This rating is based on our average peak pressure of <22.1mmHg.

Icare risk category ratings are based on average peak pressure, Therapist recommendation, case studies, industry comparison and the Waterlow scale. For some high risk clients Icare Zerotec covers are recommended. The best way to assess a client on a pressure reducing surface is to monitor pressure points with visual checks over a period of 5 days.

## CARE & CLEANING

To launder the cover in a washing machine, set the main wash to a minimum of 65°C for no less than 10 minutes. The cover can be tumble dried on low heat and the temperature must not exceed 60°C. The foam can be autoclaved at 134°C.